

Acceptable Specimen Container

Bronchoalveolar Culture



Alternate Name:		None
Performing Lab:		New Hanover
Specimen Container:		Sterile screw-cap container
Minimum Volume Required		N/A
Testing Availability	Routine: Stat:	24 hours/day No
Turnaround Time:		Preliminary reports are available after 24 hours.
Special Handling:		Specimens must be delivered to the laboratory within 30 minutes of collection.
Patient Preparation:		None
Specimen Stability:		
Reference Range:		Normal Respiratory Flora
Critical Value:		N/A
CPT Code:		87070 Culture
		87205 Gram stain
Testing Methodology:		Culture
Causes for Specimen Rejection:		Improper labeling Specimen not collected in a sterile container. Specimen container must not be contaminated on the outside. Specimens must be delivered to the laboratory within 30 minutes of
		collection.
Test includes:		Gram stain and isolation/ identification of potential pathogens, and antibiotic susceptibility testing where indicated.
Clinical Significance:		N/A