

Buffy Coat

Acceptable Tube Type



Alternate Name:		Wright Stain
Performing Lab:		New Hanover
Specimen Container:		2.5 mL lavender tube with EDTA
Minimum Volume Required:		1.0 mL
Testing Availability	Routine:	24 hours/day
	Stat:	No
Turnaround Time:		4 hours
Special Handling:		None
Patient Preparation:		None
Specimen Stability:		
Reference Range:		See CBC with Differential
Critical Value:		None
CPT Code:		85009
Testing Methodology:		Microscopic
Causes for Specimen Rejection:		Improper labeling
Other Comments:		
Clinical Significance:		N/A