

Acceptable Specimen Container

Ear Culture - Inner



Alternate Name:		Otitis media, ear culture internal
Test includes:		Isolation and identification of potential pathogens.
		Antibiotic susceptibility testing of aerobes where indicated.
		Gram stain
Performing Lab:		New Hanover
Specimen Container:		Sterile Culture Swab
Minimum Volume Required		N/A
Testing Availability	Routine: Stat:	24 hours/day No
Turnaround Time:		Preliminary reports are available after 24 hours.
Special Handling:		Specific site of specimen must be stated (Internal ear)
		Specimen must be delivered to the laboratory within 2 to 8 hours of collection. If that is not possible, the specimen must be refrigerated until delivered.
Patient Preparation:		None
Specimen Stability:		
Reference Range:		No growth
Critical Value:		N/A
CPT Code:		87070
HED Test Group:		Test Group – MLAB at NH site
HED Test Name:		Test Name – Ear Culture
Testing Methodology:		Culture

Causes for Specimen Rejection:	Improper labeling Specimen not collected in a sterile container. Specimen container must not be contaminated on the outside. Specimen not delivered to the lab within 2 – 8 hours of collection. Swabs must be received in a moist state.
Other Comments:	
Clinical Significance:	N/A