

## Acceptable Tube Type



## Free T3

Alternate Name:		FT3
Performing Lab:		New Hanover
Specimen Container:		Yellow serum separator tube
Minimum Volume Required:		1.0 mL
Testing Availability	Routine: Stat:	24 hours/day No
Turnaround Time:		Routine: 4 hours
Special Handling:		None
Patient Preparation:		None
Specimen Stability:		
Reference Range:		2.3 – 4.2 pg/mL
Critical Value:		None
CPT Code:		84481
Testing Methodology:		Chemiluminescence
Causes for Specimen Rejection:		Improper labeling
Other Comments:		
Clinical Significance:		This test is used to diagnose hyperthyroidism and to clarify thyroid status in the presence of a possible protein binding abnormality.