

Acceptable Collection Containers



**Glucose, Urine**

<b>Alternate Name:</b>		None
<b>Performing Lab:</b>		New Hanover
<b>Specimen Container:</b>		Random: Urine Specimen Cup 24 Hour: 24 Hour Urine Collection Container
<b>Minimum Volume Required:</b>		Random: 1 mL 24 Hour: Not Applicable
<b>Testing Availability</b>	<b>Routine:</b>	24 hours/day
	<b>Stat:</b>	No
<b>Turnaround Time:</b>		Routine: 4 hours
<b>Special Handling:</b>		24 Hour Collection: Requisition must state date and time collection started and date and time collection finished. Urine should be refrigerated during collection. Do not use preservatives.
<b>Patient Preparation:</b>		None
<b>Specimen Stability:</b>		
<b>Reference Range:</b>		Random: 1 – 15 mg/dL 24 Hour: < 500 mg/ 24 hr
<b>Critical Value:</b>		None
<b>CPT Code:</b>		82945
<b>Testing Methodology:</b>		Colorimetric
<b>Causes for Specimen Rejection:</b>		Improper labeling Use of preservatives
<b>Other Comments:</b>		
<b>Clinical Significance:</b>		N/A