



Rapid HIV 1 / 2 Antibody

Alternate Name:	HIV 1 / 2
Performing Lab:	New Hanover and Cape Fear
Specimen Container:	Lavender top tube
Minimum Volume Required:	5 mL
Testing Availability	Routine: 24 hours / day Stat: Yes
Turnaround Time:	Routine: 24 hours (specimens are batch tested) Stat: < 1 hour
Special Handling:	None
Patient Preparation:	None
Specimen Stability:	N/A
Reference Range:	Non-reactive
Critical Value:	N/A
CPT Code:	
Testing Methodology:	Qualitative Immunoassay
Causes for Specimen Rejection:	Improper labeling Hemolysis
Other Comments:	Positive tests are repeated in duplicate. Any repeatedly reactive patient is sent to a reference laboratory for confirmation.
Clinical Significance:	N/A