

## Acceptable Tube Type

## Intraoperative Parathyroid Hormone



Alternate Name:		iPTHOR, Intact Parathyroid Hormone
Performing Lab:		New Hanover
Specimen Container:		Lavender top tube(EDTA)
Minimum Volume Required	•	1.0 mL
Testing Availability	Routine: Stat:	24 hours/day Yes
Turnaround Time:		Routine: 4 hours Stat: < 1 hour
Special Handling:		None
Patient Preparation:		None
Specimen Stability:		2 Days at Room Temp 3 Days at 2-8° C and 6 Months -20° C
Reference Range:		15-65 pg/mL
Critical Value:		None
CPT Code:		83970
<b>EPIC Test Code</b>		LAB304178
Testing Methodology:		Chemiluminescence (Roche)
Causes for Specimen Reject	tion:	Improper labeling Hemolysis
Other Comments:		This test code is intended to be used during an intraoperative parathyroid surgery to aide the surgeon during the procedure.
Clinical Significance:		Diagnosis of parathyroid disease and other diseases of calcium homeostasis; monitoring patients undergoing renal dialysis.