

## Acceptable Tube Types

## **Intact Parathyroid Hormone**



Alternate Name: iPTH

Performing Lab: New Hanover

Specimen Container: Lavender top tube(EDTA), Yellow

serum separator tube or a red top

tube.

Minimum Volume Required: 1.0 mL

**Testing Availability** Routine: 24 hours/day

Stat: Yes

Turnaround Time: Routine: 4 hours

Stat: < 1 hour

Special Handling: None

Patient Preparation: None

**Specimen Stability:** Serum: 8 hours at Room Temp or 2-8°

C and 1 Month -20° C

EDTA Plasma: 25 hours at Room

Temp, 14 Days at 2-8° C.

Reference Range: 18.4-80.0 pg/mL

Critical Value: None

**CPT Code:** 83970

EPIC Test Code LAB108

**Testing Methodology:** Chemiluminescence (Siemens)

Causes for Specimen Rejection: Improper labeling

Hemolysis

Other Comments:

Clinical Significance: Diagnosis of parathyroid disease and

other diseases of calcium

homeostasis; monitoring patients

undergoing renal dialysis.