

New Hanover Regional Medical Center

Acceptable Tube Type





Alternate Name:		Rapid Plasma Reagin
Performing Lab:		New Hanover
Specimen Container:		7 mL red top tube
Minimum Volume Required	:	1 mL
Testing Availability	Routine:	24 hours/day
	Stat:	No
Turnaround Time:		24 hours
Special Handling:		None
Patient Preparation:		None
Specimen Stability:		
Reference Range:		Negative
		Reactive specimens will be tittered and confirmed if not previously confirmed
Critical Value:		None
CPT Code:		86592
Testing Methodology:		Charcoal Agglutination
Causes for Specimen Reject	tion:	Improper labeling Hemolyzed reactive sample
Other Comments:		
Clinical Significance:		Test for syphilis