

## Acceptable Tube Type



## Theophylline

Alternate Name:		Aminophylline
Performing Lab:		New Hanover
Specimen Container:		Red top tube
Minimum Volume Required	:	1.0 mL
Testing Availability	Routine: Stat:	24 hours/day Yes
Turnaround Time:		Routine: 4 hours Stat: < 1 hour
Special Handling:		Oral dose: Draw 4 hours after dose Slow Release: Draw 2 hours after dose
Patient Preparation:		Avoid caffeine 24 hours prior
Specimen Stability:		
Reference Range:		10 – 20 ug/mL
Critical Value:		> 20 ug/mL
CPT Code:		80198
<b>Testing Methodology:</b>		Turbidimetric – Rate
Causes for Specimen Rejection:		Improper labeling
Other Comments:		
Clinical Significance:		Theophylline is prescribed for bronchial asthma, for chronic obstructive pulmonary disease, and for newborn apnea. The drug is extensively metabolized with peak serum levels reached four hours after

oral dose.