

Acceptable Tube Types



TPN Profile

Alternate Name:	Total Parenteral Nutrition
Tests Included in Profile:	Glucose BUN Creatinine Sodium Potassium Chloride CO2 AGAP Calcium Phosphorus Magnesium Triglyceride Prealbumin Albumin
Performing Lab:	New Hanover and Cape Fear
Specimen Container:	Yellow, Red, or Green top tube or Microtainer
Minimum Volume Required:	1 mL
Testing Availability	Routine: 24 hours/day Stat: Yes
Turnaround Time:	Routine: 4 hours Stat: 1 hour
Special Handling:	None
Patient Preparation:	None
Specimen Stability:	N/A
Reference Range:	See individual tests
Critical Value:	See individual tests
CPT Code:	See individual tests
Testing Methodology:	See individual tests
Causes for Specimen Rejection:	Improper labeling

Hemolysis

Other Comments:

Clinical Significance:

See individual tests